

Compact DR plus

Motorized unit with batteries



CONTACT US

INTERMEDICAL S.r.l.
IMD Group

24050 Grassobbio (Bergamo) Italy
Via E. Fermi, 26
ph. 0039 035 6594811
fax 0039 035 6594899
e-mail info@inter-med.it
www.inter-med.it

Data may be altered with no prior warning.

doc. 11394 - 28/05/2013



IMPROVING PERFORMANCES

Mobile Unit equipped with digital images acquisition system in radiography mode with detector at solid status.

**DIGITAL VERSION
WITH FLAT PANEL**



ANALOGIC VERSION



COMPACT DR, suits a wide range of clinical applications, intended to perform diagnostic and radiological examinations of:

- Intensive care
- First Aid
- Orthopaedics
- Neonatology
- Sports Medicine
- Operating Theatre

The acquired images are saved in the system memory, displayed on a high resolution monitor and possibly sent to a diagnostic work station or to a general archive through Dicom connection.

The unit grants the automatic calibration of the acquisition parameters, process and image display according to both the selected anatomic part and the patient size.

**Wide vertical excursion
(569mm ÷ 2069 mm)**

1. Cable reel
2. Electric and mechanic arm blocking mechanism
3. Safety key and DVD recorder
4. Anticollision bumper



1.

2.

3.

4.

