


User-friendly interface

Softwares available:

- **DIM** (RAM or Flash memories)
Option: USB (Dicom compatibility)
- **AFG** and **DFG** memories
(HD with 350.000 images DVD, USB)
Modular configurations even after sales
Options: Cine, Measures, Road Mapping, DSA, DICOM

Washable keyboard
with trackball




CONTACT US

INTERMEDICAL S.r.l.
IMD Group

24050 Grassobbio (Bergamo) Italy
Via E. Fermi, 26
ph. 0039 035 6594811
fax 0039 035 6594899
e-mail info@inter-med.it
www.inter-med.it

Data may be altered with no prior warning.


doc. 10289 - April 2013




Radius
C-arm




IMPROVING PERFORMANCES



200 mm horizontal C-arm run,
with manual brake for locking



Orbital movement 125°
Free space: 840 mm
SID:1040 mm



270° on each side arm rotation,
with manual brake for locking.



12° on each side C-arm swivelling
with manual brake for locking.



*A reliable, tough and enduring
C-arm with the maximum
affordability; intuitive software
for different configurations with
user-friendly interface.*

- *Functionality*
- *Software modularity*
- *Image acquisition – processing
– management – storage*
- *Lower dose*
- *Better features*